

A Spectrum of Products with Just What You Need

EUTECH
INSTRUMENTS
Technology Made Easy ...



pH | ORP | Ion | Conductivity | TDS | Dissolved Oxygen | Temperature

Eutech 700 Bench Meters—Simple, Affordable Solutions for Water Testing

You shouldn't have to pay for features you don't use. The Eutech 700 Bench Meter Series offers five meters with the features you need for simple, affordable measurements of pH/ORP/Ion, Conductivity/TDS, Dissolved Oxygen and pH/ORP/Conductivity/TDS—all with °C/°F Temperature.

- Largest display of any bench meter in it's class—more than twice the size of our popular 510 series.
- Simple, stress-free meter operation.
- High accuracy, reliable meters feature auto temperature compensation, non-volatile memory, auto-buffer recognition, splash-proof keys and more!
- Multi-position electrode arm and quick reference guide included.



Chemopharm Sdn Bhd (25504-W)

20, Jalan SS2/66, 47300 Petaling Jaya, Selangor, Malaysia

T +603-7872 6000 | F +603-7874 7958

enquiry@chemopharm.com.my | www.chemopharm.com.my

For more information, please contact:

Eva Yeong @ 012-278 1365

Scott Liew @ 012-396 8196