

*nova*  
Blood Glucose  
Monitoring System  
**Max™**  
Quick. Simple. Smart.



- ✓ **No Coding Needed**
- ✓ **World's Smallest Blood Sample**
- ✓ **Fast 5-Second Results**
- ✓ **Alternate Site Testing**

### Nova Max™ Monitor

The Nova Max Monitor offers accuracy and precision in blood glucose testing, with more convenience than ever before.



**Quick**

Fast Test Time – 5 seconds

**Simple**

Small Blood Sample – 0.3 microliters Alternate Site Testing – Decrease pain associated with testing

**Smart**

No Coding Needed – Accurate and precise tests 400 Test Memory – Review past test results with ease

### Nova Max Link™ Monitor

The Nova Max Link Monitor offers accuracy and precision for pump patients.



**No Coding Needed**

Miscoded monitors can lead to Insulin dosing errors which could be as high as 50% 1.2

**Proven Technology**

Second generation monitor using proven technology

**Communicates to Medtronic Insulin Pump**

Nova Max link automatically sends new blood glucose test results to the Paradigm Insulin pump for use by the Bolus Wizard Calculator.

**Hassle-Free Upgrade**

Paradigm Link users can upgrade without a new prescription for test strips.

For more information about the Nova Max Link, please visit the website <http://www.novamaxlink.com>

### Software

By using Nova Diabetes Software's colorful and easy-to-read charts, you and your health care provider will be able to see the "patterns" when your blood glucose values are in and out of target.



To download free software and purchase cable, please visit the website <http://www.novacares.com/software/software.html>

### Accessories

By using Nova Diabetes Software's colorful and easy-to-read charts, you and your health care provider will be able to see the "patterns" when your blood glucose values are in and out of target.



**Nova Sureflex Lancing Device**



**Nova Sureflex Lancets**



**Nova Max Control Solution**